

# Viewing Digital Images

Digital Images in pathology lab reports are viewed in an application called Techcyte Fusion AP, which is accessible via MayoACCESS.

To view Digital Images:


1. If you are not already logged in to [MayoClinicLabs.com](https://www.mayocliniclabs.com) and with an active MayoACCESS session in progress, click the **MayoACCESS** hyperlink in the website banner, then log in.



For detailed login instructions, please see [Getting started](#).


2. Access UniView:

- a. From the **Results** menu, select **Pathology Lab Reports** (right).



The **Pathology Lab Reports** page appears, showing the **Reports** tab by default (left). The displayed list of available reports defaults to the last 2 weeks.

- b. To change search parameters, click the **Search Criteria** tab at top left, complete the appropriate fields, and, at top right, click the **Search** button.



c. Select the desired order row in your search results to highlight it.

Pathology Lab Reports

TEST,IHC ID: C7 8-IHC Female 01/01/2001 24Y Prim. Phys: SSN:

Reports Search Criteria

UT Unsolicited Tests RI Result Images

\* Report times for Mayo performed tests are CST/CDT

Patient	ID	Order #	Accession #	Collected	Reported*	Ab	Img	Status
TEST,IHC	C7 8-IHC	M181388081	M181388081	02/14/23 01:00	02/14/23 14:28			Final

Report Result List View Supplemental Data

Mark All Reports As Read


d. To view images, perform one of the four options below:

- **Option I:** In the SmartLinks Bar or on the SmartMenu, click the **Report** SmartLink.


The first screenshot shows a SmartLinks Bar with the following options: Report, Result List, Fax Report, View Supplemental Data, and Mark All Reports As Read. The 'Report' link is highlighted with a red box.

The second screenshot shows a 'Pathology Lab Reports SmartMenu' with the following options: Export Grid, Grid Report, Fax Report, Mark Report As Read, Report, Mark All Reports As Read, Refresh Grid, Result List, and View Supplemental Data. The 'Report' link is highlighted with a red box.


- **Option II:** In the list of reports, double-click the selected order row.  
Either Option I or II loads the report in the MayoACCESS **Report** window.





In the report, find the test you want to view images and/or data for, and click the link below **Supplemental data is available at:**




- **Option III:** In the selected order row, in the **Img** column, click the paper icon (right).




- **Option IV:** In the SmartLinks Bar or on the SmartMenu, click the **View Supplemental Data** SmartLink.






Either Option III or IV brings up the **Supplemental Data for Order** dialog box (left).



In the row of the test you want to view, click the link in the **Supplemental Data URL** column (left).

Techcyte Fusion AP opens in a new browser tab. The **Slides** icon (screenshot below, upper left corner) is selected by default and displays the Pathology Image Window.



For instructions on using Techcyte Fusion AP, see the [application user guide](#).