

## Viewing Digital Images


Digital Images in pathology lab reports are viewed in an application called Techcyte Fusion AP, which is accessible via MayoLINK

To view Digital Images:

1. Access Fusion AP by one of the following methods:

- **Option 1:**

a. Click the link in the **Test Results Available** email notification.



b. Log in to MayoLINK. For instructions, see [Logging in to MayoLINK](#).

The MayoLINK **Orders** page opens in a new browser tab. Only the order linked in the email notification appears in the Orders List.

- **Option 2:**

- If you are not already logged in to [MayoClinicLabs.com](https://www.mayocliniclabs.com) with an active MayoLINK session in progress, click the **MayoLINK** hyperlink in the website banner, then log in.

For detailed login instructions, see [Logging in to MayoLINK](#).

- On the **Select Account** page in MayoLINK, choose the number of the account for which you want to view digital images.
- Find the desired order by using the MayoLINK **Search** or **Advanced Search** function or scrolling through the Order List.

For detailed instructions on using MayoLINK search functionalities, see [Searching for an order](#).

2. View images by clicking one of the two **Status** column links described and illustrated below:

The screenshot shows the MayoLINK Orders interface. At the top, there are navigation tabs: Orders, Batches, Patients, Preferences, Metrics, Audit Log, and Help. Below this is a search bar with a search button and a search filter dropdown. The main content area displays a table of orders. The first row is highlighted in yellow and contains the following information:

Order Number	Patient Information	Tests Ordered	Status
MAYOLINK-092823-1 70 28 Sep 2023 15:29	CONSULTANTS, RENAL K UAT-RENAL1	• Renal Pathology Consultation, Wet Tissue + Request More Tests	Display Report PDF JPG TIF PNG TXT Result Details Revised

A red box highlights the Status column for the first order. A red arrow points from this box to a larger, detailed view of the Status column options. This detailed view shows the following options:

- Display Report (with a document icon) - Option II
- PDF
- JPG
- TIF
- PNG
- TXT
- Result Details (with a monitor icon) - Option I
- Final


- **Option I:** Click the monitor icon.
- **Option II:** Click the **Display Report** link to open the report PDF. Then, on the report, click the link under **Digital images available at:** (below).

The screenshot shows a Mayo Clinic Laboratories report page. At the top, there is the Mayo Clinic Laboratories logo and the phone number 1-800-533-1710. Below this is a patient information table:

Patient ID DIGITALIMAGES	Patient Name SECTRA, PATHOLOGY S	Birth Date 2012-12-18	Sex U	Age 10
Order Number MQ00051828	Client Order Number 091323-TR-KR	Ordering Physician [REDACTED]	Report Notes	
Account Information C70	Collected 13 Sep 2023 14:17			

Below the patient information table, there is a section titled **BOB-1 IHC, Tech Only** with a sub-section **Performed**. The text states: "The positive and negative control is verified as showing appropriate immunoreactivity." Below this, there is a section titled **Digital images available at:** with a link: <https://qa.orders.mayocliniclabs.com/en/slides/> Sectra UniView Case: KR-23-214. A red arrow points to this link. Below the link, there is a section titled **Renal Pathology** with a sub-section **Interpretation**. The text states: "Received: 13 Sep 2023 14:47" and "Reported: 13 Sep 2023 15:56". Below this, there is a section titled **FINAL DIAGNOSIS** with a sub-section **COMMENT**. The text states: "Case Number KR-23-214".

Techcyte Fusion AP opens in a new browser tab. The **Slides** icon (screenshot below, upper left corner) is selected by default and displays the Pathology Image Window.



For instructions on using Techcyte Fusion AP, see the application user guide.