

## Viewing Digital Images


Digital Images in pathology lab reports are viewed in an application called Techcyte Fusion AP, which is accessible via MayoLINK

To view Digital Images:

1. Access Fusion AP by one of the following methods:

- **Option 1:**

a. Click the link in the **Test Results Available** email notification.



b. Log in to MayoLINK. For instructions, see [Logging in to MayoLINK](#).

The MayoLINK **Orders** page opens in a new browser tab. Only the order linked in the email notification appears in the Orders List.

- **Option 2:**

- If you are not already logged in to [MayoClinicLabs.com](https://www.mayocliniclabs.com) with an active MayoLINK session in progress, click the **MayoLINK** hyperlink in the website banner, then log in.


For detailed login instructions, see [Logging in to MayoLINK](#).

- On the **Select Account** page in MayoLINK, choose the number of the account for which you want to view digital images.
- Find the desired order by using the MayoLINK **Search** or **Advanced Search** function or scrolling through the Order List.

For detailed instructions on using MayoLINK search functionalities, see [Searching for an order](#).


2. View images by clicking one of the two **Status** column links described and illustrated below:

The screenshot shows the MayoLINK 'Orders' page. At the top, there are navigation tabs: Orders, Batches, Patients, Preferences, Metrics, Audit Log, and Help. Below this is a search bar with a 'Search' button and a 'Download Reports' button. A filter bar shows 'All', 'Incomplete (43)', 'Pending Batch', 'Sent to Lab', 'Reported (50+)', and 'Results Processed'. A 'Download list' link is visible. The main table has columns for Order Number, Patient Information, Tests Ordered, and Status. The Status column for the first order is highlighted with a red box, showing options: 'Display Report' (with a PDF icon), 'PDF', 'JPG', 'TIF', 'PNG', 'TXT', 'Result Details' (with a monitor icon), and 'Revised'. A red arrow points from this box to a larger, detailed view of the 'Status' column options below. In this detailed view, 'Display Report' is labeled 'Option II' and 'Result Details' is labeled 'Option I'.

- **Option I:** Click the monitor icon. 
- **Option II:** Click the **Display Report** link to open the report PDF. Then, on the report, click the link under **Digital images available at:** (below).

The screenshot shows a report from Mayo Clinic Laboratories. The header includes the logo and 'MAYO CLINIC LABORATORIES' with the phone number '1-800-533-1710'. The report details include Patient ID (DIGITALIMAGES), Patient Name (SECTRA, PATHOLOGY S), Birth Date (2012-12-18), Sex (U), and Age (10). The test name is 'BOB-1 IHC, Tech Only'. The report states 'Performed' and 'The positive and negative control is verified as showing appropriate immunoreactivity.' Under 'Digital images available at:', there is a link: <https://qa.orders.mayocliniclabs.com/en/slides/...> with 'Sectra' and 'UniView Case: KR-23-214'. The report was received on 13 Sep 2023 14:47 and reported on 13 Sep 2023 15:56. The case number is KR-23-214. The report includes sections for 'Renal Pathology', 'Interpretation', 'FINAL DIAGNOSIS', 'COMMENT', and 'MICROSCOPIC DESCRIPTION'.

Techcyte Fusion AP opens in a new browser tab. The **Slides** icon (screenshot below, upper left corner) is selected by default and displays the Pathology Image Window.



For instructions on using Techcyte Fusion AP, see the [application user guide](#).