

If you select the **Digital Image** option when creating your rule, the email notification (below, left) will contain a link to a report that itself contains a link to the images. The recipient clicks the link in the report (below, top right), and Techcyte Fusion AP opens, displaying the relevant images (below, bottom right).

The image is a composite of three parts illustrating a workflow:

- Top Left (Email):** An email notification from Mayo Clinic Laboratories. The subject is "INT Test Results Available for 7 1 001*****017q". The body contains a link to a report: "001*****017q".
- Top Right (Report Snippet):** A snippet from a report titled "BOB-1 IHC, Tech Only". It includes patient information: Patient ID SA00000179, Patient Name Miss UNSOLICITED-ENVTWO, TECHCYTE R, Birth Date 2002-04-15, Sex F, Age 23. A red box highlights a link: "Supplemental data is available at: https://mci-int-portal.mayo.edu/accounts/C7234577/orders/3388388_order-4nd5-078a-3881484b5c98/tech/RER17/images".
- Bottom (Techcyte Fusion AP):** A screenshot of the Techcyte Fusion AP interface. The URL is "https://techcyte.com/cases/...". The case ID is "TR-26-80000". The interface shows a large image of a tissue slide with brown spots, a sidebar with "Slides" (A: 1-PR-DfDx uTV8i, B: 1-BOB1 uTAs), and a right sidebar with "Case details" and "Annotations".